

SCALING

SCALING



MODEL G1

- · Removal of supra gingival tartar.
- Removal of tartar, narrow interdental space.
- · Removal of marginal gingival tartar.









ORDER CODE Z291101

Z217101

Single rooted / Multi rooted

Cross Section Diagram

EMS® MODEL G2-E

MODEL G2

• Removal of supra gingival tartar.

23.4

ORDER CODE Z291102

Z217102

Cross Section Diagram

NSK Satelec® MODEL G3 EMS® MODEL G3-E

MODEL G3

• Removal of supra gingival tartar.

20.0

Z217103 ORDER CODE Z291103

Cross Section Diagram

MODEL G4

- Removal of supra gingival tartar.
- Removal of marginal gingival tartar.







Z217104 ORDER CODE Z291104 MODEL G5

- · Removal of supra gingival tartar.
- Removal of tartar, narrow interdental space.
- · Removal of marginal gingival tartar.
- Removal of gingival biofilm under natural tooth.











- · Removal of supra gingival tartar.
- Removal of tartar, narrow interdental space.
- · Removal of marginal gingival tartar.
- Removal of gingival biofilm under natural tooth.









Z217106 ORDER CODE Z291106

Cross Section Diagram















MODEL G8

Removal of supra gingival tartar.

Removal of tartar, narrow interdental space.

Removal of marginal gingival tartar.

22.7 Cross Section Diagram

MODEL G9

Removal of supra gingival tartar.

• Removal of tartar, narrow interdental space.

• Removal of marginal gingival tartar.

Removal of gingival biofilm under natural tooth.



MODEL G11

Removal of supra gingival tartar.



-612 22.9 Cross Section Diagram



NSK Satelec® MODEL G8

Z217108

Z217109

NSK Satelec® MODEL G11

Z217411

NSK Satelec® MODEL G11

MODEL G12

• Removal of supra gingival tartar.

Z217412

MODEL G13

· Removal of supra gingival tartar.

NSK Satelec® MODEL G11

Z217413

Cross Section Diagram











* Cross section diagram is the image of a point 1 mm from the end. The illustrated tip diagram is four times bigger than the actual size.

Refer to the power guide on pages 55-56.